

# Biomanufacturing with the *i*CapTag™: a Novel Scalable and Tagless Protein Purification Platform




PROTEIN PURIFICATION TECHNOLOGY  
CONTACT US & COLLABORATE

Izabela Gierach, Ph.D., MBA  
CEO, Co-Founder  
Protein Capture Science

David W. Wood, Ph.D.  
CSO, Co-Founder  
Protein Capture Science  
Professor at the Dept. of Chemical and Biomolecular Engineering  
The Ohio State University

## Background

The purification of recombinant proteins at laboratory and manufacturing scales represents a substantial difficulty in the development of new therapeutic proteins. At laboratory scales, researchers rely on affinity tags to rapidly purify diverse new targets. These approaches cannot be used at manufacturing scales however, due to the potential immunogenicity of the tags and the difficulty of removing the tags at manufacturing scales. The result is a divided world of protein purification, where tag-based methods dominate in the laboratory, but non-tag methods must later be developed for manufacturing. The disconnect between these methods can significantly delay drug development and approval, and in some cases can stop promising therapeutics from ever being commercialized. The end result is untreated diseases, high drug costs, and slow market entry for new lifesaving medicines.



**Figure 1.** Image of Npu intein, a main component of patented *i*CapTag technology, which is able to self-assemble under specific pH-buffer conditions also self-cleave allowing on target protein to be seamlessly purified. The two pieces of Npu are shown in **Figure 2** as a red box cartoon, that can bind together, "disassemble" during stripping step and reassemble again during target protein binding, where target protein exist as a complex and is expressed together with Npu C-terminal. On another hand cleavage of protein can be triggered by pH change.

## Definitions

**Inteins** (Intervening proteins) sequence is located between two exteins, N- and C-terminals, in nature. Interestingly, inteins naturally can splice or cleave, which led to development of intein-based purification technology here described. Moreover, when rationally engineered inteins can self-cleave and self-assemble in predictable way. Inteins-based tag technology here presented is enabled through a C-terminal cleavage reaction induced simply by pH shift of the buffer. Note, that intein cleavage rate is pH and temperature dependent.

## Acknowledgement

Thank you graduate students working on development and testing of early prototypes of protein purification technologies, including Changhua Shi, Joe Taris (Streptokinase and G-CSF work here presented), Jackelyn Gallardi, Brian Marshall, Tzu-Chiang Han, and Yamin Fan. Special thanks goes to funding agencies, including NSF, NIH, DARPA, Ohio Third Frontier, OSU Fund, and industry partners.

## Method


The *i*CapTag technology is a single platform for protein purification and is designed to be used in small and large scale. Therefore, steps here described are reminiscent to those seen in ordinary R&D lab developing new biopharmaceuticals and in manufacturing (see **Figure 2**).

### FLEXIBILITY OF EXPRESSION SYSTEMS

The *i*CapTag can be used with commonly used expression systems such as mammalian and bacterial host-cells to develop biopharmaceuticals in R&D labs as well as by using cell-free expression systems. Based on internal data, the tag here described does not interfere with the secretion of expressed target protein. At the same time, it was noted that this technology is ideal for the proteins that can be well expressed and fold.

### NEUTRALITY TO OTHER BUFFER ADDITIVES


While the pH is controlled the *i*CapTag technology was found neutral to salts, sugars and other buffer additives including protease inhibitors used for decrease of proteolytic degradation.



**Figure 2** Steps of protein purification using *i*CapTag technology that bridge R&D and manufacturing.

### CLEAVAGE RATE CONTROLLABILITY & PREDICTABILITY

The cleavage rate (k) of the target protein of interest is simply controlled by shift of pH from e.g., pH 8.5 to pH 6.2. Based on internal research, the rate of cleavage and pH sensitivity reminds mostly the same at pH 8.5 independently of first "-1" position of amino acids (a.a.) located at the beginning of Npu N-terminal. However, pH sensitivity and cleavage rate dramatically increase at pH 6.2 and depend on addition of "-1" type of a.a. At the same time, "+1" position of a.a. located at the beginning of target protein (added or naturally occurring) was found to be critical for controllability of cleavage rate (smaller & charged a.a. led to slower k, aromatic & larger a.a. led to faster k).



### ABILITY OF RESIN REGENERATION

The *i*CapTag technology can be regenerated by stripping buffer containing 150 mM phosphoric acid and 500 mM sodium chloride (pH  $\leq 2.0$ ).

### STERILIZATION BUFFERS

The *i*CapTag can be used with caustic buffers without a loss of resin activity.

### ICAPTAG APPLICATIONS



The *i*CapTag can be used for novel and well-known proteins, peptides, fragments of proteins and biosimilars. This technology is also ideal for high-throughput studies, and it can be easily implemented within existing R&D labs and CMO/CRO set-up. This technology is also scalable and can be used as a single platform from R&D to manufacturing.

### TAGLESS PROTEIN VIA ONE PLATFORM

It was shown that the beauty of the intein-based technology is that the cleavage of the target protein leads to precise and tagless release of protein. It was also shown that in some cases, the activity of the tested proteins such as G-CSF (blue line) was higher than standard purchased from PeproTech (see red line shown in **Figure 3**).

## Results

Examples of protein purifications of known biopharmaceuticals using intein-based technology are shown below for R&D usage only.



## References

1. Boffert M, Boffert G, Delychew N, Wood DW, Wu W. Inventors. Reversible Polypeptide Intein. assignee: Genetic system and self-cleaving intein derived therefrom, bioprocess and protein purification employing same, and methods for determining critical gene-sizable amino acid residues for varying intein activity. US patent 8,933,362. August 24, 2020.
2. The Protein Data Bank. H. A. Berman, J. Westbrook, Z. Feng, G. Gilliland, T. N. Weiss, M. Weissig, J. D. Weissig, J. D. Weissig. Nucleic Acids Research. 28: 235-242.
3. Wu WY, Miller KD, Cookson M, Wood DW. Inteins-mediated one-step purification of Escherichia coli secreted human antibody fragments. Protein expression and purification. 2011;76(2):120-126.
4. Berman TO, Cookson M, Wood DW. Ligand-independent control and self-cleaving intein as a tool for high-throughput protein purification. Protein expression and purification. 2012;92:100-104.
5. Shi C, Han TC, Wood DW. Purification of Microbially Expressed Recombinant Proteins via a Dual ELP-Split Intein System. Methods Mol Biol. 2017;1485:1-25.
6. Wu W, Miller KD, Shi C, Han TC, Wood DW. Intein-mediated control of protein expression and purification. Protein expression and purification. 2011;76(2):120-126.
7. Lohs A, Fan Y, Stoppel SD, Rahn M, Wood DW. Intein as a tool for tagless and tagless protein purification. Journal of Chemical Technology and Biotechnology. 2017;97(1):167-180.
8. Fan Y, Miller KD, Stoppel SD, Han TC, Wood DW. Column-Free Purification Method for Recombinant Proteins Using Self-Cleaving Tagging. Protein. 2018;11(5):464-474.
9. Fan Y, Miller KD, Stoppel SD, Han TC, Wood DW. Column-Free Purification Method for Recombinant Proteins Using Self-Cleaving Tagging. Protein. 2018;11(5):464-474.
10. Izabela Gierach and David Wood. Self-Cleaving Tags Based on Split Inteins Increased Reliability Enabling Higher Throughput Applications. American Pharmaceutical Review. 2019; 22(1): 86-89.